




How to use Braun ThermoScan® 5



Remove the thermometer from the protective case.



Press the Power button.



Attach a new lens filter.



Fit the probe snugly into the ear canal.




Press and release the Start button.



Temperature displays.

• Always replace disposable lens filters to maintain accuracy and hygiene.



Product	NEW IRT 6510 IRT 6030	
Patented pre-warmed tip	✓	✓
Positioning system	✓	✓
Ready beep	✓	✓
Memory storage capacity	9	1
Easy-to-read LCD display	✓	✓
Protective case with lens filter storage	✓	✓
1.5V AA LR6 2x	✓	✓
Disposable lens filters 21x	✓	✓
Color	Blue	Green

Braun Healthcare Malaysia
www.braunhealthcare.com

These products are manufactured under license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany. ThermoScan® and ExacTemp® are trademarks of Helen of Troy Limited and/or its affiliates.

© 2018 All rights reserved.

Distributor:
Pharmaforte (Malaysia) Sdn. Bhd.
2, Jalan PJU 3/49, Sunway Damansara,
47810 Petaling Jaya,
Selangor Darul Ehsan,
Malaysia.
Tel: 603-7803 6998



Watch product video

Aug 2018



BRAUN



The ear thermometer pediatricians recommend the most¹.

ThermoScan® 5
Infrared ear thermometer



IRT 6510 / IRT 6030



BRAUN

**Do you know how to
take my temperature as
precisely as a doctor ?
My thermometer does.**

When it comes to accurately taking a child's temperature, doctors trust the precise Braun ThermoScan[®] infrared ear thermometer.

ThermoScan[®] provides professional accuracy thanks to its patented, pre-warmed tip. Disposable Braun lens filters help to eliminate the risk of cross-contamination. And with its large screen and ergonomic design, ThermoScan[®] is simple to use. Discover why the Braun ThermoScan[®] is the No.1 brand among pediatricians¹. Doctors know best.

Features



More accurate

Exclusive pre-warmed tip technology is designed to minimize the cooling effect of the tip inside the ear canal, providing accurate measurements time after time.



Unique positioning system

ExacTemp[™] technology uses a light and beep to confirm proper positioning.



Night light feature

Helps take a temperature even when it's dark.



Memory function

-IRT 6510: Recalls last 9 temperature readings.
-IRT 6030: Recalls last temperature reading.



Safe and hygienic

Braun disposable lens filter to eliminate the risk of cross-contamination and ensure accurate readings.

Clinically proven more accurate

than rectal or temple³.

How it works

Braun ThermoScan[®] measures the infrared heat generated by the eardrum and surrounding tissue. Because the eardrum shares blood supply with the temperature control centre in the brain, measurement with Braun ThermoScan[®] accurately reflects changes in core body temperature.

Guaranteed accuracy for all ages

With its patented pre-warmed soft tip and guidance system, the Braun ThermoScan[®] infrared ear thermometer is suitable for babies ears from birth onwards.

Sources:

1. Based on a US survey reporting that of the surveyed pediatricians (1) who recommend a brand of thermometers, Braun is the most recommended; and (2) who use a branded thermometer, Braun is the most frequently used.
2. Based on research conducted by Nielsen, Jan-Feb 2017, in five Asia Pacific markets [n=1465 Mom of kids 0-12 years old in AU, HK, KO, SG, TW]
3. Measurements in detecting temperature changes. Nimah, Marianne M., et. al. "Infrared tympanic thermometry in comparison with other temperature measurement techniques in febrile children." *Pediatr Crit Care Med*. 7.1 (2006): 48-55.

What is fever?

If your body temperature is consistently above your normal temperature range,* you are said to have fever. Fever is a defense mechanism the body uses to fight infections. However, fever alone isn't an illness, it's a symptom.

* Practice taking temperatures on yourself to determine the normal temperature range.

Would the tip on the thermometer affect accuracy reading?

Yes. The tip on other ear thermometers can cool the area where the reading is taken, leading to a higher risk of inaccurate readings. The patented Braun pre-warmed tip is designed to minimize the cooling effect of the tip inside the ear canal, providing accurate measurements time after time.

How can parents reduce the risk of infection when their children have a fever?

Infants and children can suffer from fever many times throughout the year, and parents usually take their child's temperature several times during each bout of fever, depending on how unwell the child feels.

Taking measurements with Braun disposable lens filters enable hygienic use without having to clean or disinfect the sensor and virtually eliminate cross-contamination. What's more, fitting a new Braun lens filter each time ensures that dirt or ear wax don't build up and clog the filter membrane, which could lead to inaccurate temperature readings.